

Radiation Safety. Amplified.

Tyo-Cycle II

BERRA

Cryo-Cycle[™] II

Hybrid Cryostat







Nuclear

Homeland Security

& Defense

Labs and Industrial and Education

Manufacturing

KEY FEATURES

- LN₂ redundancy
- Non-CFC/non-flammable refrigerant
- Low power demand
- Same footprint as standard LN₂ Dewar

Healthcare

- Long life Pulse-Tube cooler Lifetime (L5) >75000 hours of continuous operation
- Remote read-out and control
- Low vibration/low electrical noise
- Available in dipstick and integral configurations
- 2-year full warranty + pro-rated • warranty on the cooler

BENEFITS

- Low operating cost
- Higher up-time
- Field installable (dipstick version)
- Quiet •
- No compromise on detector specifications

DESCRIPTION

The CANBERRA Cryo-Cycle™ II is a unique offering in the field of cryogenically cooled radiation detectors. The Cryo-Cycle II is described as a "hybrid" cryostat because it combines the advantages of electric cooling with the reliability of liquid nitrogen. The long-life Pulse-Tube cryocooler, condenses the boil-off N₂ gas back into the 25 liter Dewar. This unique capability provides the convenience of operating a detector for 12 to 18 months before LN₂ needs to be added, but at the same time keeps the detector cold in case of power failure. With the Cryo-Cycle II the LN₂ supply keeps the detector cold for up to one week without power. There is no interruption of cooling. There is no downtime due to partial warm-up as long as LN₂ level is maintained. There is no risk of detector failure because of temperature cycling. LN₂ lost during power outages may be replenished at any time.

Cryo-Cycle

CANBERRA

The Cryo-Cycle II comes with a number of improvements allowing us to answer our customers' requirements even better.

Key Improvements

- Reduced cooler vibrations
- Improved LN₂ level sensor probe
- Single front panel
- RS-232 and USB interfaces
- Alarm and autofill relay outputs
- A Graphical User Interface (GUI)

The audible noise has been reduced to less than 60 dB(A), measured at 1 m distance, making the Cryo-Cycle II well suited for application in guiet laboratory environments.

The new LN₂ level sensor probe provides better accuracy. The measured LN₂ level is displayed through a continuous LED indicator scale on the front panel, allowing to better schedule periodic refills.

CANBERRA

All controls, connectors and indicators are integrated in a single front panel for easy access. The Cryo-Cycle II contains an auto-ranging power supply at 100-240 V and 50-60 Hz. The front panel is equipped with a DB9-M, a RS-232 and a USB connector. The DB9-M is a relay output for system status alarms and autofill functions. The serial connectors are used to connect to a PC. A dedicated GUI allows remote control and status monitoring.



Cryo-Cycle II front panel layout

The GUI is available through a Windows[®]-based software application, provided with the Cryo-Cycle II. This application needs to be installed on a PC connected to the Cryo-Cycle II through the USB or RS-232 serial ports. Minimum operating system requirements are Windows XP (SP3) or Windows 7.

Cryo-Cycle Control Panel			
File Commands Help			
LN Level		Status	
	Cooler	OFF	
	Cooler Power	0 Watts	
100	Pressure	0 psi	
80	Heater	OFF	
50	Autofill	disabled	
40	Autom	disabled	
30		Alarms	
	No	Alarm Conditions	
Not connected		User mode	Log disabled

Screenshot of the Cryo-Cycle II control panel application

The application can be operated in User mode or in Supervisor mode. The User mode only allows status and parameter monitoring, while the Supervisor mode, which can be password protected, allows access to the available commands. The displayed parameters and status indications can be continuously logged to a user-selectable .txt-file, saved on the PC's hard drive. Also each Cryo-Cycle II can be given a system name, allowing easy identification when multiple systems are monitored through the same PC. This name will be displayed in the application's title bar.

The Cryo-Cycle II is designed to accommodate both dipstick and integral configurations. Dipstick versions can be installed in the field, while integral versions must be assembled at the factory.

Due to the improved microphonics performance of the Cryo-Cycle II, when it is sold with a new CANBERRA detector, there will be NO degradation of the detector's resolution performance as stated on the detector's specification sheet. If the dipstick version is installed on older CANBERRA detectors some degradation of resolution performance may occur, depending on the age and configuration of the detector. CANBERRA guarantees no resolution degradation at energies above 500 keV and a maximum of 10% between 100 and 500 keV. Performance is not guaranteed below 100 keV. For detectors not manufactured by CANBERRA, resolution performance cannot be guaranteed.

The highly reliable and efficient Pulse-Tube cooler (lifetime of >75000 hours) used in the Cryo-Cycle II contains a CFC free and non-flammable gas. The cooler is hermetically sealed, so no gas-refill is required. The compressor contains no oil or lubricant, so no contamination of the refrigerant occurs and no periodic filter/dryer exchange is required. This makes the Cryo-Cycle II virtually maintenance free. The nominal power consumption is very low (250 W), with a maximum of 450 W in transient operation. The Cryo-Cycle II is designed to operate between 10 °C and 35 °C.

CANBERRA's confidence in the Cryo-Cycle II is demonstrated by the two year full warranty on the complete system (detector included when purchased together) and an additional pro-rated warranty on the cooler. If the cooler fails after the second year, it will be repaired or replaced at 40, 60 or 80% of the list price in year three to five respectively. This pro-rated warranty applies to parts only.

SPECIFICATIONS

PERFORMANCE

- CANBERRA guarantees detector performance as warranted by detector model with cooler in operation (on new detectors purchased with Cryo-Cycle II).
- LN₂ loss rate <3 liters/day typically (with cryocooler OFF).
- MAINTENANCE Cleaning as required to keep air flow unrestricted.
- LEVEL INDICATORS Linear LED scale on front panel.

CONNECTORS

- USB 2.0 Remote control and status read-out.
- RS-232 Remote control and status read-out.
- DB9-M Relay output.

COOLING

• Forced air (internal fans).

POWER REQUIREMENTS

- 100–240 V ac, 50–60 Hz, 690 VA max. (auto ranging power supply).
- FUSE (2) T 5 A 250 V (193-240 V ac Operation).
 FUSE (1) T 10 A 250 V (100-130 V ac Operation).
- NOMINAL POWER CONSUMPTION -250 W.

PHYSICAL

COLD HEAD (EXCLUDING DETECTOR CHAMBER)

- DIMENSIONS 43.2 cm (17 in.) diameter x 61.0 cm (24 in.) high.
- WEIGHT 28.2 kg (62 lb) empty, without detector.
- DEWAR-CAPACITY 25 liters.

ENVIRONMENTAL

- OPERATING TEMPERATURE +10 to +35 °C (50 to 95 °F) on standard models and configurations.
- OPERATING HUMIDITY RANGE: 20% to 80% relative non-condensing.
- Meets the environmental conditions specified by EN 61010, Installation Category I, Pollution Degree 2.

SOFTWARE

- SYSTEM REQUIREMENTS Windows XP (SP3) or Windows 7 (32-bit).
- .NET framework 3.5 (will be installed if not present, requires internet connection).

AVAILABLE DETECTOR MODELS

 Cryo-Cycle II can be ordered with all standard GC-, GX-, GR-, BE-, and GSW-detector models (see applicable detector specification sheets for details).

ORDERING INFORMATION

Model	Description
CCII-VD	Cryo-Cycle II for model 7500SL or 7500
CCII-HD	Cryo-Cycle II for model 7600SL or 7600
CCII-VI-SL	Cryo-Cycle II vertical integral Slimline
CCII-VI-F	Cryo-Cycle II vertical integral Flanged
CCII-HI-SL	Cryo-Cycle II horizontal integral Slimline
CCII-HI-F	Cryo-Cycle II horizontal integral Flanged
CCII-HI-U	Cryo-Cycle II horizontal U-Style

CRYO-CYCLE II WITH VERTICAL DIPSTICK CRYOSTAT (CCII-VD)



Cryo-Cycle II with 7500SL cryostat



End cap dimensions depend on detector size. The tables below show the typical surface area or efficiency range vs. end cap diameter. End cap lengths are also greater for larger detectors. Consult the factory if end cap size is critical in your application.

LEGe/BEGe, Nom. Area (mm ²)	End Cap Diameter, mm [in.]
=<2000	76 [3.0]
2800	83 [3.25]
3800	89 [3.50]
5000	102 [4.0]
6500	114 [4.50]

Coax Rel. Efficiency (%)	End Cap Diameter, mm [in.]
=<40	76 [3.0]
40-50	83 [3.25]
50-70	89 [3.50]
70-100	95 [3.75]
100-120	102 [4.0]
120-150	108 [4.25]
150	114 [4.50]

Cryo-Cycle II with 7500 cryostat







CRYO-CYCLE II WITH HORIZONTAL DIPSTICK CRYOSTAT (CCII-HD)

Cryo-Cycle II with 7600SL cryostat





Cryo-Cycle II with 7600 cryostat



CRYO-CYCLE II WITH HORIZONTAL INTERAL CRYOSTATS



Cryo-Cycle II Horizontal Integral U-Style (CCII-HI-U)













CRYO-CYCLE II WITH VERTICAL INTEGRAL CRYOSTATS











Cryo-Cycle is a trademark and/or registered trademark of Mirion Technologies, Inc. and/or its affiliates in the United States and/or other countries.



©2017 Mirion Technologies (Canberra), Inc. All rights reserved.

Copyright ©2017 Mirion Technologies, Inc. or its affiliates. All rights reserved. Mirion, the Mirion logo, and other trade names of Mirion products listed herein are registered trademarks or trademarks of Mirion Technologies, Inc. or its affiliates in the United States and other countries. Third party trademarks mentioned are the property of their respective owners.

CANBERRA

