

# PRESSURE LOSSES IN BENDS AND FITTINGS HB100D



**1** study

## **Features**

- · A self standing demonstration unit
- Unit contains various pipe types and of differing diameters
- Connections of differing types
- · Uses valves to control flow

### **Description**

The unit demonstrates the pressure losses that can occur in different types of pipes and plumbing arrangements. The unit comprises of 6 different horizontal sections, different diameter, roughness and material pipes, along with long and short radius bends, parallel sections and constrictions. The unit also has an expandable top section for looking at different available options or students own designs.

Isolating valves are fitted to each of the sections so that

each component can be examined individually. A hand held digital manometer fitted with quick release couplings allows quick and accurate measurement of the pressure drop across any of the pipes/fittings.

In conjunction with the HB100K the unit can be used to examine various different pipe flow theorems covering topics such as Friction factors, Reynolds numbers, flow velocity's and laminar boundary layers for turbulent pipe flow in each of the different sections.

#### **Related laws**

- · Reynold's number
- Friction Factor
- · Bernoulli's equation
- Fluid Velocity
- Laminar Boundary layer
- Plumbing



- Hydraulics
- · Mechanical engineering
- · Chemical engineering

#### Learning capabilities

- Determination of the pressure drop across a variety of pipe sections, materials and sizes
- Determination of Reynolds number in a variety of pipe sizes
- Determination of friction factor for a variety of pipe sections materials and sizes.
- Calculation of theoretical pressure loss in a pipe using Bernoulli's equation – comparison with practical measured pressure drop.
- · Determination of fluid velocity in pipes and fittings
- · Determination of laminar boundary layer thickness

## **Technical Specification**

- 90° Elbow
- · Parallel pipes
- 45° swept tee
- · Long elbows
- Short elbows
- Sudden constriction
- 1/2" dia copper pipe
- 20mm dia PVC pipe
- 32mm dia rough bore PVC pipe.
- 32mm dia smooth bore PVC pipe.
- 10 valves.
- 22 pressure tapping points.
- Expandable top section with 3 separating modules.
- Flexible inlet and outlet allowing unit to be rotated 90° for space saving.
- Detachable feet to allow wall mounting.

#### **Essential Ancillaries**

- HB100/230
- or
- HB100/115
- HB100L\_230

#### What's in the Box?

- Experimental Module
- Drain Hose Assembly
- Water Inlet Straight Pipe Assembly
- Removable Frame Feet Assembly with Locking Castors
- Digital Manometer and Hoses
- Water Inlet 90° Elbow Assembly

## Weights & Dimensions

Weight: 30 kgLength: 2200mmWidth: 150mm

#### **Essential Services**

· Height: 1000mm

- HB100/230
- or
- HB100/115
- HB100K

### **Ordering information**

To order this product, please call PA Hilton quoting the following code: HB100D

All brand and/or product names are trademarks of their respective owners. Specifications and external appearance are subject to change without notice. The colour of the actual product may vary from the colour shown in the brochure.

Copyright © 2018 P.A. Hilton Limited. All rights reserved. This technical leaflet, its contents and/or layout may not be modified and/or adapted, copied in part or in whole and/or incorporated into other works without the prior written permission of P. A.

**COUNTRY OF ORIGIN - UK WARRANTY PERIOD - 2 YEARS** 

Hilton Limited. Hi-Tech Education is a registered trade mark of P. A. Hilton Limited.