



Waste treatment by absorption

Référence : CAB/1000



Solution decontamination by a reaction with a gas

Continuous operation

GENERAL SPECIFICATIONS

- 2 tank 30 L.
- Column, analyzer with sheath pH, plate with sampling, ozone inlet plate with sintered.
DP measurement.
- Dosing pump.
- Anti-pulsing pot.
- Filter regulator with isolation valve.
- Hydraulic guard.

Instrumentation

- Flowmeter.
- Gas flowmeter .
- Manometer.

Dim : 105 x 55 x 215 cm – 100 kg
SS tubular framework 40 x 40mm

Categories: Educational Engineering Environment Waste treatment

Reference: CAB/1000

DESCRIPTION

A method of gas-liquid absorption at different water flow rates

Example: study of ozone disinfection



Unités livrées avec un manuel pédagogique
et dossier technique



Possibilité d'adapter les unités à vos
besoins



Mise en service, formation des
formateurs